




## Verifying the HUMAC2015 Servo Settings

### 1. Introduction


- 1.1. The HUMAC System is compatible with a number of different isokinetic testing machines and servo amplifiers. This document describes the procedure to set the HUMAC program for the correct machine and servo amplifier.
- 1.2. After installing the HUMAC on a new computer you must start and exit the HUMAC program one time before running this procedure.

### 2. 64-bit HUMAC Interface - NORM

- 2.1. From the **Windows Task Bar**, select **Start, All Programs, HUMAC2015, Utilities, Diagnostics64**.
- 2.2. From the **Startup Screen**, enter **CSMIDIAG** as the password and click the **Proceed with Caution** button.
- 2.3. The HUMAC will display the **HUMAC64 Diagnostics** form. The message *"Initializing HUMAC Interface Board. Please wait..."* will be displayed at the bottom of the page. During this time the HUMAC is loading the **PLD** and **DSP**, and reading the **MACHINE**, **SERVO**, and **MOTOR** settings.



- 2.4. When the procedure is complete, the **MACHINE SETTINGS** area will display the **MACHINE**, **SERVO**, **MOTOR** and **REV** values.



- 2.5. Click the **EXIT** button to save the settings.


**3. 32-bit HUMAC Interface.**

- 3.1. From the **Windows Task Bar**, select **Start, All Programs, HUMAC2015, Utilities, Diagnostics**.
- 3.2. From the **Startup Screen**, enter **CSMIDIAG** as the password and click the **Proceed with Caution** button. The HUMAC will display the **USB-DSP-PLD monitor** form.
- 3.3. The **USB-DSP-PLD monitor** form will display the current Machine, Servo, Motor and Rev settings.
- 3.4. Follow the appropriate section for the NORM, 6000, Biodex System2 and System3 machines.


Machine	Section
NORM	4
SYSTEM2	5
SYSTEM3	6

**4. NORM**

- 4.1. Set the Machine to NORM.
- 4.2. Leave the Rev set to 0.




- 4.3. In the **HUMAC Tools** area, click the **Load PLD** button. The HUMAC will load the PLD. This takes about 10 seconds.




- 4.4. In the **Servo Programming** area, click the **Get Servo** button. This procedure will take about 10 seconds.



- 4.5. The HUMAC will display the current Machine, Servo and Motor settings. In this case the **NORM, ULTRA3000** and **Getty**.




- 4.6. If the system displays **DM-30**, **DM-75** or **Ultra3000**, the procedure is completed. Click the **Exit** button to end the procedure.




- 4.7. If the system displays “**UNKOWN**”, repeat step 4.4 through 4.6 one more time. If the system still displays “**UNKOWN**”, proceed to the section titled “**Manually Setting the NORM Servo Amplifier**”.

4.8. Manually Setting the NORM Servo Amplifier


- 4.8.1. In the **HUMAC Tools** area, click the **Load DSP** button. This will take about 5 seconds.








- 4.8.2. Click the **Start DSP** button.  
 4.8.3. Click the **HUMAC Tool** button.  
 4.8.4. Remove all input adapters from the NORM.  
 4.8.5. Open both ROM stops.  
 4.8.6. Position the input adapter with the arm at 6:00 (Straight Down).  
 4.8.7. In the **Supervisory Commands** area, click the **Reset** button.




- 4.8.8. Wait 5 seconds and click the **Self Test** button.  
 4.8.9. Rotate the dynamometer arm three complete revolutions in the clockwise direction, ending at the 6:00 position.  
 4.8.10. Record the **Position** value displayed (17308 in this case).



- 4.8.11. Click the “**X**” in the Upper-Right-Hand corner to close the **HUMAC Tool**.  
 4.8.12. On the right-hand side, you will see three pull-down lists. Set them as follows.

Recorded Position Value	<i>Action Required</i>
<p style="text-align: center;">Near 17300</p> 	<p>System is set correctly.</p>
<p style="text-align: center;">Near 16000</p> 	<p>System is using an ULTRA3000/Getty but is incorrectly set to.DM30 or DM 75 and Rockwell.</p> <ol style="list-style-type: none"> <li>1. Select the following options Ultra 3000 and Getty.</li> </ol>  <ol style="list-style-type: none"> <li>2. Repeat steps 4.8.3 thru 4.8.10 to confirm the position value is now near 17300.</li> </ol>
<p style="text-align: center;">Near 18888</p> 	<p>System is using a DM30 or DM 75 but is incorrectly set to ULTRA3000/Getty.</p> <ol style="list-style-type: none"> <li>1. Open the Servo Housing and determine if the NORM has a DM30 (4" wide) or DM75 (5-1/2" wide) installed.</li> <li>2. Select the correct servo and set the Motor to Rockwell.</li> </ol>  <ol style="list-style-type: none"> <li>3. Repeat steps 4.8.3 thru 4.8.10 to confirm the position value is now near 17300.</li> </ol>

4.8.13. If the value near 17300 is displayed, click the Exit button to end the procedure.




**5. Biodex System2**

- 5.1. Set the **Machine** to **System2**.
- 5.2. The **Servo** and **Motor** will be set to **IPC** and **QMC**.
- 5.3. Set the **Rev** as follows:

<b>Double Chair (Rev 0)</b>	<b>Single Chair (Rev 1)</b>
-----------------------------	-----------------------------

Machine: SYSTEM2	Machine: SYSTEM2
Servo: IPC	Servo: IPC
Motor: QMC	Motor: QMC
Rev: 0	Rev: 1

5.4. The procedure is completed. Click the **Exit** button to end the procedure.



**6. Biodex System3**

- 6.1. Set the **Machine** to **System3**.
- 6.2. The **Servo** and **Motor** will be set to **IPC** and **QMC**.
- 6.3. Set the **Rev** as follows:

Red dot on Knee/Hip Adapter is in-line with adapter post (Rev 0)	Red dot on Knee/Hip Adapter is rotated 90 degrees from adapter post (Rev 1)
Machine: SYSTEM3	Machine: SYSTEM3
Servo: IPC	Servo: IPC
Motor: QMC	Motor: QMC
Rev: 0	Rev: 1

6.4. The procedure is completed. Click the **Exit** button to end the procedure.

