


| | | | |
|---|---------------------------------------|---|------------------------|
|  | Computer Sports Medicine, Inc. | | |
| | Title: | Determining Your NORM Dynamometer Revision | |
| | Document #: | 300051 | |
| | Rev: | A | Effective Date: |

CYBEX produced two types of NORM dynamometers. The first design used two black plastic connectors at the back of the dynamometer. The second design used three metal connectors at the back of the dynamometer. To determine the type of NORM you have, look at the wires leading out of the dynamometer. The first type has two round plastic connectors. The second type has three round metal connectors. Please indicate the type of NORM you have.

Facility Name: _____

Phone: _____

_____ I have a two black plastic connector dynamometer.



Figure 1 - Two Connector NORM

_____ I have a three metal connector dynamometer.



Figure 2 - Three Connector NORM

Please complete and FAX this sheet to CSMI at 781-297-2039.